

IN STOCK!

CE MEDICAL FACE MASK - TYPE IIR

BFE 98 % - SPLASH RESISTANCE $\geq 16,0$




Table 1 — Performance requirements for medical face masks

Test	Type I ^a	Type II	Type IIR
Bacterial filtration efficiency (BFE), (%)	≥ 95	≥ 98	≥ 98
Differential pressure (Pa/cm ²)	< 40	< 40	< 60
Splash resistance pressure (kPa)	Not required	Not required	$\geq 16,0$
Microbial cleanliness (cfu/g)	≤ 30	≤ 30	≤ 30

^a Type I medical face masks should only be used for patients and other persons to reduce the risk of spread of infections particularly in epidemic or pandemic situations. Type I masks are not intended for use by healthcare professionals in an operating room or in other medical settings with similar requirements.

DON'T SAVE ON QUALITY

					
	Chirurgisch masker	FFP2 FFP3 zonder ventiel	FFP2 FFP3 met ventiel	Stoffen masker	Sjaal
Beschermt uzelf	WEINIG	JA	JA	WEINIG	WEINIG
Beschermt anderen	JA	JA	NEE	WEINIG	WEINIG
Hospitaal-gebruik	JA	JA	JA	NEE	NEE

